

CONFIDENTIAL B

The MediaTek logo consists of the word "MEDIATEK" in white, uppercase, sans-serif font, centered within an orange parallelogram shape that is wider at the top and bottom and tapers in the middle.

**MEDIATEK**

# MediaTek IoT SmartDevice App Introduction

**2017.06.10**

# Outline

- Overview
- SmartDevice Features
  - User Scenario
  - BTNotify
  - Connection
  - Features (Notification Push/FOTA/FMP/BAS/EPO)
  - About
- Health Features
  - Overview
  - Health Feature


# Overview

# SmartDevice App (1/3)


- **MediaTek SmartDevice App** is an Android application used for MediaTek IoT device (based on **MT2523/MT2533** chip).
- It doesn't only apply to MediaTek IoT Bluetooth device (MT2523/MT2533), but also for Smart Health Device (MT2523 + MT2511, watch\_ref\_design project).

# SmartDevice App (2/3)

## ■ App GUI



(Connected with IoT Device,  
disable Health Feature)



(Connected with 2511 health device,  
enable Health Feature)


# SmartDevice App (3/3)

- The APK includes below features:
  - Scan nearby BT Devices
  - Switch Mode (SPP/GATT)
  - Connection (Based on MTK BTNotify transport protocol)
  - Notification Push
  - FOTA (Firmware Over-The-Air by BTNotify transport)
  - EPO Download (only applies for MTK IoT GNSS Project)
  - Find Me & BAS (only enable in GATT mode)
  - MT2511 Health Feature (only enable for health device)
    - Heart Rate Measure
    - Blood Pressure Measure
    - Blood Pressure Calibration & Personal Mode
    - Heart Rate Multi-connect & Comparison

# SmartDevice Features


# SmartDevice Main UI

- Main Activity







# User Scenario



MainActivity  
Scan nearby BLE Device



DeviceScanActivity  
Each scan only 60 seconds




MainActivity  
Connect successfully

# BTNotify (1/2)

- BTNotify
  - BTNotify: MTK BT Transport Protocol
  - Handshake: APK must handshake to confirm that the remote device is MTK Smart Device (based on BTNotify) after connect BT device successfully.
  - Two Mode
    - SPP (Based on BT SPP Profile)
    - GATT (i.e. DOGP Mode, Based on MTK defined BLE GATT Profile - DOGP - Data over GATT Profile)


# BTNotify (2/2)



# Connection (1/3)

## ■ Switch Mode


- Support SPP and GATT mode.
- SPP and GATT cannot coexist.
- APK will connect only BTNotify SPP profile in SPP mode, but will not connect BT A2DP/HFP profiles.
- APK will connect BLE profiles (BTNotify DOGP/FMP/BAS) in GATT mode.
- APK can switch mode in disconnected state.
- APK cannot switch to GATT mode if SP don't support BLE.
- APK will keep the last mode when SP or APK reboot.



**Current mode is SPP, APP will switch to GATT mode after click “Switch Mode”.**

# Connection (2/3)

- Scan
  - Show EDR/DUAL devices in SPP mode.
  - Show bonded device in SPP mode directly.
  - Show LE/DUAL devices in GATT mode.



# Connection (3/3)



- Auto Reconnect
  - GATT/SPP auto reconnect the last GATT/SPP connection address, GATT/SPP will auto connect after reboot BT, SP, Remote Device.
  - GATT/SPP doesn't auto connect after user disconnect GATT/SPP connection by click operation. (i.e. Code call disconnect API)

# SmartDevice Features

- SmartDevice main features:
  - Notification Push
  - FOTA
  - FMP
  - BAS
  - EPO

# Notification Push (1/3)

- Enable:
  - If you don't activate Notification Push, APK will ask you to activate it in System Setting UI.






**Activate Notification Push**



# Notification Push (2/3)

- Notification App List:
  - You can disable the APP that if you don't want to sync its notification to device.



# Notification Push (3/3)

- SmartDevice will listen new arrived info/notification and push them to remote device, such as:
  - **Notification**: Normal android application notification. User can block some app's notification in APK. So when receive new notification of this app, the APK will don't push it to remote device.
  - **Message**: If receive new SMS, APK will push it to remote device.
  - **Missed call**: If miss a call, APK will notify to remote device.
  - **Low battery warning**: If SP changes to low battery status, APK will notify to remote device.

# FOTA (1/3)

## ■ FOTA


- Firmware Over-The-Air , send Firmware data to remote device for Firmware Upgrade.

## ■ Workflow

- Get current version from remote device via BTNotify.
- Choose firmware upgrade file from local file system.
- Transfer firmware data via BTNotify.
- Remote device starts to upgrade after transfer completely.
- Feedback the upgrade result to APK after remote device reboot and BTNotify auto reconnect successfully.

# FOTA (2/3)



- FOTA Flow



# FOTA (3/3)

## ■ FOTA UI

- FOTA menu will be available after connected successfully.
- FOTA could send full device bin, or GNSS bin(only for GNSS project) to remote device.



Only demo,  
Not real  
firmware file




# FMP (1/3)

- Find Me Profile (FMP)

- When a button is pressed on one device to cause an alerting signal on a peer device.


- Find Phone

- Press button on watch, the phone will alert out




- Find watch

- Press button on Phone, the watch will alert out



# FMP (2/3)

## ■ FMP UI




SP APK initiates “Find Me” alert/stop to find device.

(Only enable in MT2533 Device)

# FMP (3/3)

- FMP UI




Remote Device initiates “Find Me” alert to find Smart Phone.  
(Enable in MT2523/MT2533 related project)




# BAS (1/2)

- The **Battery Service** (BAS) exposes the state of a battery within a device
  - SP Read battery level from device when connected.




- When device battery level change, it should notify phone.



# BAS (2/2)

- BAS UI
  - SmartDevice APP working as BAS Client.



(Only enable in MT2533 Device)



# EPO

## ■ EPO

- EPO: A file be used to assistant location in GNSS.
- The EPO file will be download from network server and send to device once request and SP has been connected to the internet.
- Only enable in MT2523 GNSS project.

# About (1/2)


- About GUI



SmartDevice APK send a notification to testing notification push feature.

# About (2/2)

## ■ About GUI



SmartDevice APK send a demo command by using DemoController.

# Health Features

# Overview (1/3)

## ■ MT2511

- Bio-AFE support PPG (photoplethysmography) / ECG (electrocardiography) signal.
- MT2523/2533 can communicate with MT2511 bio-sensor via I2C and SPI.

## ■ IoT Device SW

- Use IoT SDK “*project\mt2523\_watch\apps\watch\_ref\_design*”
- Feature
  - HR (PPG + Motion)
  - Blood Pressure (PPG + ECG)
  - Touchgfx GUI

# Overview (2/3)


## ■ APP Health Feature

- Add record – Transfer personal profile to device
- HR Measure – measure heart rate
- BP measure
  - General Mode – Measure personal blood pressure
  - Personal Mode – Measure bp with saved personal model, which created from calibration process.
- Multiple connect – Add reference device for heart rate comparison



# Overview (3/3)

## ■ APP GUI




# Add Record

- Add Record
  - Fill in the personal profiles and transfer to device, which will be used when measuring HR/BP.
  - Note: This is mandatory option when measuring BP.

The screenshot shows a mobile application interface for adding a new record. The title bar at the top is brown and contains a back arrow, the text 'BP Measure', a document icon, and a plus sign. The main content area is a white form titled 'Add Record' in orange. The form contains the following fields: 'User ID' with the value 'Jack' (highlighted in blue), 'Age' with the value '29', 'Gender' with 'Male' selected (indicated by an orange circle) and 'Female' as an option (indicated by a grey circle), 'Height (cm)' with the value '173', and 'Weight (kg)' with the value '68'. At the bottom of the form are two buttons: 'CANCEL' in grey and 'OK' in orange. Below the form, at the very bottom of the screen, is a brown button labeled 'START'.


# HR Measure (1/2)

- Present bio-sensor data
  - Entry or Turn on heart rate measurement on your device, you will see PPG chart, HR BPM and confidence level on the screen.
  - Confidence level
    - 0 ~ 1: low confidence level
    - 2 ~ 3: high confidence level
  - If you would like to edit the personal profile, click on the “+” (upper right of the screen).




# HR Measure (2/2)

- Store raw data
  - Tap on the “START” button on the bottom of the screen will start a recording session.
  - Tap on the same button (which should be showing “STOP” now) again will stop the recording session and save the result to the local storage on the smart phone. (/sdcard/CatchLog).




# BP Measure – General Mode (1/3)

- Present bio-sensor data
  - Turn on blood pressure measurement on your device and wait for a few seconds, you can see both PPG and ECG chart on the measurement area.
- Store raw data
  - Tap “START” button will start a recording session.
  - If you have a golden device, use the second button to the upper right to enter the golden measurement result as a record for debugging usage.




# BP Measure – General Mode (2/3)

- Present blood pressure result
  - Show systolic BP, diastolic BP, and HR BPM.
  - Show PWTT chart.




# BP Measure – General Mode (3/3)

- SNR
  - If user store raw data (“Start” recording), SNR10 for PPG and SNR40 for ECG would be computed and be shown on screen.
  - User can evaluate the signal quality with the SNR values.



# BP Measure – Personal Mode (1/6)

- In personal mode, user can either creates a new personal model or chooses a personal model which had been stored in APK.
- User will execute a calibration process when he/she decided to create a personal model.





# BP Measure – Personal Mode (2/6)

- There are 4 steps in calibration process,
  1. Create a personal model
  2. Input Golden
  3. Calibration
  4. Input Golden again
- User can create his own personal model through these steps.
- The **first** step is simple, just like general mode, user has to input profiles such as birth year, gender, height, and weight.

Calibration

**1. Create a personal model**  
In order to create your personal model, please input your profile.

User ID  
mediatekB

Birth Year  
1980

Gender  
☒ Male ☐ Female

Height (cm)  
180

Weight (kg)  
50

CANCEL OK

Start


# BP Measure – Personal Mode (3/6)

- During calibration process, user need to use a sphygmomanometer as a golden device to measure his/her current blood pressure.
- In the **second** step, input the measure result from the golden device.

The screenshot shows the 'Calibration' screen of the 'mediatekB' app. The screen has a dark grey background with a grid pattern. At the top, there's a brown header bar with the word 'Calibration' and a user icon. Below the header, the text 'mediatekB' is visible. A white modal box is centered on the screen with the title '2. Input Golden' and the instruction 'Please enter the measure result from sphygmomanometer.' Inside the modal, there are two input fields: 'Systolic BP' and 'Diastolic BP'. At the bottom right of the modal are 'CANCEL' and 'OK' buttons. Below the modal, the text 'REC time 00:00:00' is displayed. At the very bottom of the screen is a large, rounded 'Start' button.


# BP Measure – Personal Mode (4/6)

- In **third** step, user would switch features to select BP measure or calibrate mode.
- When user selects it, the measure of personal model will start and measuring personal model for 25 seconds.
- During the measurement process, EKG and PPG raw data will be transferred to SP through BT and the waveforms will be shown in APK screen.




# BP Measure – Personal Mode (5/6)

- **Finally**, after the personal calibration parameter had been measured and transferred to SP, user will be asked to input the measured result from golden device again.
- When user press OK, all the personal model data will be stored in APK.
- The calibration process is completed and user can measure blood pressure with their own model.




# BP Measure – Personal Mode (6/6)

- When user chooses a personal model to measure blood pressure, his personal model would be set to the device.
- User can measure his blood pressure based on the personal model which he had created before through calibration process.
- Just like the process in general mode.



# Multiple Connect (1/4)

- Multiple connect feature
  - Select the reference device in the list and tap the “CONNECT” button.
  - Make sure the device is discoverable and support SIG BLE HR Profile.




# Multiple Connect (2/4)

- Multiple connect feature
  - Heart rates from both devices will now show in the measure screen.
  - Tap the “START” button to start a recording session.


MT2511 Health Device.

Reference HR Device.



# Multiple Connect (3/4)

- Multiple connect feature
  - The AP5 and AP10 percentage during the recording period will be calculated and shown on the screen.





# Multiple Connect (4/4)

- Multiple connect feature

- Formula:

- Diff Count (total number about our hr bpm is difference from golden bpm be measured from reference device)
    - AP5 Count (total number of our hr bpm differ 5~10 from golden bpm)
    - AP10 Count (total number of our hr bpm differ 10 or higher from golden bpm)
    - AP5:  $1 - (\text{AP5 Count}) / (\text{Diff Count})$
    - AP10:  $1 - ((\text{AP5 Count}) / (\text{Diff Count})) * ((\text{AP10 Count}) / (\text{Diff Count}))$

# Contact US

(If you have any questions, comments, or suggestions, please contact us by MTK ACS, or send mail to [SmartDevice\\_App@mediatek.com](mailto:SmartDevice_App@mediatek.com))



*everyday genius*